

Mobile X-ray System IMAX 1010 V



IMAX 1010V is a lightweight, portable device that can be mounted on a mobile stand and used as a mobile ward X-ray unit.

The touchscreen software enables easy configuration of anatomical programs (APR)



FireCR 20 is a digital radiography device that uses phosphor plate technology. It is a cost-effective solution offering high resolution and throughput.

Technical Specifications

Generator power	8 kW
Tube voltage	40 – 125 kV
Focal spot size	0.6 mm / 2.8 mm
Sampling pitch	Selectable: 200 µm or 100 µm
Cassette sizes	14" × 17" and 14" × 10"

Configuration

Item name	Quantity
-----------	----------

Mobile X-ray System IMAX 1010 V:

- | | |
|--|---|
| <ul style="list-style-type: none">• Stand• Generator• X-ray tube• VacuDAP compact dosimeter | 1 |
|--|---|

Workstation with FireCR 20 Digital Radiography Device:

- | | |
|--|---|
| <ul style="list-style-type: none">• Digital radiography device• Two cassettes: 14" × 17" and 14" × 10"• 23" monitor• Software | 1 |
|--|---|

e-mail: medical.trest@gmail.com

+37062812800