



Rescue Life Defibrillator

Rescue Life is a new leading product developed through the improvement of Progetti defibrillators. With over 25 years of experience, Rescue Life represents an advanced solution for emergency medicine. Its ergonomic and extremely durable design makes it the first professional defibrillator fully manufactured in Italy.

The automatic self-test mode allows checking the proper functioning of the device, the battery status, and the charge-discharge system, with real-time data management and report printing capabilities.

Designed to meet the needs of modern emergency medicine, Rescue Life features a large 7-inch color TFT LCD display, a powerful battery providing up to 5 hours of continuous monitoring, and multiple configuration options.

Rescue Life Defibrillator



- Electrical protection class: Class I
- Battery operating time: More than 5 hours
- Battery capacity: 150 discharges at 230 J
- Memory: 4 GB internal memory with up to 300 hours of recording
- Patient age group: Adults / Children
- Water protection: IPX4
- Waveform: Biphasic truncated exponential
- Patient impedance: From 25 to 200 Ohm
- Display: 7-inch TFT with backlight
- Display resolution: 800 × 480 pixels
- Output energy: From 1 to 230 J
- Electrodes: Adult and pediatric defibrillation via detachable or attachable reusable external paddle contact plates
- Operating mode: Manual / Semi-automatic
- Synchronization: Synchronized to the “R” peak
- ECG gain: 2.5 mm/mV, 5 mm/mV, 10 mm/mV, 20 mm/mV, 40 mm/mV, Auto
- Lead selection: 5 leads (10 optional)
- Heart rate (HR): Measurement range 20–300 bpm

Scope of Supply – Rescue Life Defibrillator

- Defibrillator main unit – 1 pc
- Reusable adult and pediatric electrodes – 1 set
- Built-in printer – 1 pc
- Lithium battery – 1 pc
- 5-lead ECG cable – 1 pc
- AED / Advisory module with:
 - Voice prompts
 - 2 GB internal memory
 - 1 cable
 - 1 set of disposable adult pads
 - – 1 set
- User manual – 1 pc

Rescue Life Defibrillator



- Electrical protection class: Class I
- Battery operating time: More than 5 hours
- Battery capacity: 150 discharges at 230 J
- Memory: 4 GB internal memory with up to 300 hours of recording
- Patient age group: Adults / Children
- Water protection: IPX4
- Waveform: Biphasic truncated exponential
- Patient impedance: From 25 to 200 Ohm
- Display: 7-inch TFT with backlight
- Display resolution: 800 × 480 pixels
- Output energy: From 1 to 230 J
- Electrodes: Adult and pediatric defibrillation via detachable or attachable reusable external paddle contact plates
- Operating mode: Manual / Semi-automatic
- Synchronization: Synchronized to the “R” peak
- ECG gain: 2.5 mm/mV, 5 mm/mV, 10 mm/mV, 20 mm/mV, 40 mm/mV, Auto
- Lead selection: 5 leads (10 optional)
- Heart rate (HR): Measurement range 20–300 bpm

Scope of Supply – Rescue Life Defibrillator

- Defibrillator main unit – 1 pc
- Reusable adult and pediatric electrodes – 1 set
- Built-in printer – 1 pc
- Lithium battery – 1 pc
- 5-lead ECG cable – 1 pc
- AED / Advisory module with:
 - Voice prompts
 - 2 GB internal memory
 - 1 cable
 - 1 set of disposable adult pads
 - – 1 set
- User manual – 1 pc

Rescue Life Defibrillator



Installation & Commissioning

- Installation, start-up, and calibration are performed by the company's engineers

Staff Training

- On-site training on how to operate the system
- Online training provided by the manufacturer from Italy
- Training on working with the software is also included

Warranty Period

- 12 months from the date of commissioning
- An additional post-warranty service contract can be signed if required