

COMMERCIAL OFFER

Digital Electrocardiograph E65



Features:

- 12-channel electrocardiograph
- 8-inch color LCD display
- Power consumption: 100 VA
- Continuous operation up to 5 hours
- Charging time: no more than 6 hours
- Memory for 800 ECG data groups
- Heart rate measurement: 30–300 bpm
- Real-time arrhythmia detection
- Real-time waveform freeze (“freeze frame”) function
- Rechargeable lithium-ion battery
- Operating modes: automatic, manual, and rhythm analysis
- External memory support: flash drive, SD card
- Automatic background data storage
- ECG report conversion to PDF format
- Selectable filter system depending on the type of signal interference
- Defibrillation protection
- Paper speed: 5 mm/s, 10 mm/s, 25 mm/s, 50 mm/s
- Dimensions (W × H × D): 360 mm × 276 mm × 130 mm
- Weight: 4.2 kg



The electrocardiograph is intended for use in medical facilities.

Equipment комплектation:

- Main device unit
- ECG cable
- Suction cup electrodes – 6 pcs
- Clip electrodes – 4 pcs
- Thermal paper – 1 roll (144 sheets, 210 × 140 mm, 20 m)
- Power supply cable

Manufacturer:

Guangdong Biolight Meditech Co., Ltd., PRC (China)