

## **Operating Table 300 B**

Manufacturer: Atese Sağlık Hizmetleri Ltd. Şti. (Turkey)

**Designed for surgical, orthopedic, proctological, gynecological, otolaryngological (ENT), urological, and other similar procedures.**

- Provides maximum access to the patient during surgical intervention.



**Maximum access to the patient during surgical procedures.**



- Convenient positioning for the patient, surgeon, nursing staff, and anesthesiologist.
- Simple and intuitive operation.
- High patient load capacity, regardless of the patient's weight and height.
- Eco-friendly and sterile design.

### Table Positions – Technical Specifications (EN)

- Head section: 55° up · 90° down, removable
- Lateral tilt: 25° right · 25° left for side tilt positioning
- Trendelenburg position: 30°
- Reverse Trendelenburg: 25°
- Back section: 60° up / -40° down ( $\pm 15^\circ$ )
- Leg section: 30° up · 90° down · 90° with lateral sliding, detachable
- Type: Electromechanical (controlled)
- Longitudinal table movement for use with C-arm X-ray system: 300 mm



### Technical Specifications

- Table top: Radiolucent
- Minimum height: 710 mm
- Maximum height: 1010 mm
- Table top dimensions: 600 mm × 2150 mm
- Head section: 225 mm × 355 mm
- Back section: 530 mm × 625 mm
- Seat base section: 530 mm × 530 mm
- Right leg section: 255 mm × 710 mm
- Left leg section: 255 mm × 710 mm
- Maximum load capacity: 400 kg
- Power supply: 220 V, 50 Hz
- Built-in battery: Yes
- Battery : 45 minutes



# MEDICAL TREST

## **Delivery Time**

Within 35 working days across the entire territory of Ukraine from the date of contract signing.

## **Installation & Commissioning**

Installation, start-up, and calibration are performed by the company's certified engineers.

## **Staff Training**

On-site training on system operation is provided at the installation location.

## **Warranty Period**

12 months from the date of commissioning.

## **Documentation**

The system undergoes full state certification and is supplied with a user manual.