

MICRONOMIC® III (Entermed B.V., The Netherlands)



MICRONOMIC® III is the embodiment of ergonomics, functionality, and modern design for the otolaryngologist. It is specially designed for intensive daily use in clinical environments, including hospitals, outpatient departments, private practices, and medical centers.



Key Advantages

- **Compact design** — suitable for limited-space environments.
- **High mobility** — easy to move thanks to four lockable castors.
- **Intuitive operation** — automatic sensors for quick function activation.
- **Superior hygiene** — integrated containers for instruments and medical waste.
- **Flexible configuration** — modular system with the option for future expansion.
- **European quality** — CE certified, manufactured in the Netherlands.



System Configuration

Cabinet construction:

Laminated multilayer structure finished in Signal White (RAL 9003)

Work surface:

Lacquered with textured finish (Signal White, RAL 9003)

Front panels:

Metal with textured coating White Aluminium (RAL 9006)

Storage drawers:

Two drawers for instrument storage with internal dividers

Small aluminum tray:

For mounting on an articulated holder

Large aluminum tray:

Integrated into the tabletop

Containers:

With openings for used instruments and waste

Power supply:

1 × 230V socket for connection of additional equipment

Configuration:

Right-sided functional arrangement

Mobility:

Mounted on 4 lockable castors

Integrated Systems

Aspiration System (without secretion container)

Includes:

- Aspiration handpiece with tubing
- Reservoir with overflow protection

Features:

- Silent operation and maintenance-free
- Power: 35 L/min
- Vacuum: -0.80 bar


LuminENT III Lighting System (LED)

- LED light source with fiber-optic output
- Compatible with STORZ cable
- Automatic ON/OFF via sensor in the handpiece

Thanks to its compact dimensions,

high-quality materials, and easy access to all functions, MICRONOMIC® III significantly optimizes the workflow of the ENT specialist while ensuring maximum safety and comfort for both the physician and the patient.

MICRONOMIC® III is a mobile diagnostic and treatment unit with a modular configuration concept, ideally suited for limited working spaces. Its smooth mobility on castors and well-organized layout of functional elements guarantee intuitive operation and efficient performance of medical procedures.



Technical Specifications

Power Supply:

230 V ($\pm 10\%$), 50 Hz

Dimensions (W × D × H):

87 × 90 × 74 cm (minimum 81 × 66 cm)

Weight (with integrated systems):

100 kg

Color:

RAL 9003 (White), RAL 9006 (Aluminium)

ENT Chair TENT 4000 (TISCOMED, Turkey)

A modern fully automated chair designed for examinations and minor procedures in otorhinolaryngology. It combines comfort, ergonomic design, and high functionality to ensure optimal working conditions for the physician and maximum comfort for the patient.



Technical Specifications

Control system:

Fully automated system

Rotation capability:

360°

Position adjustment:

Up/down, forward/backward

Footrest:

Movable

Upholstery:

Soft padding made of special molded polyurethane foam, upholstered in high-quality synthetic leather (antibacterial, fire-resistant, wear-resistant artificial leather)

Chair control:

Touch control panel or foot switch

Return-to-start function:

Automatic return to the initial position with one button

Frame construction:

All metal parts of the body are coated using an electrostatic painting method

Total Offer Price

Item Name	Qty
MICRONOMICÂ® III	1
ENT Chair TENT 4000	1

Delivery, Installation & Warranty

Delivery time:

Within 40 business days across the entire territory of Ukraine from the date of contract signing.

Installation & commissioning:

Installation, start-up, and calibration are performed by certified company engineers.

Staff training:

On-site system operation training is provided at the installation location.

Warranty period:

24 months from the date of commissioning.

Documentation:

The system undergoes full state certification and is supplied with complete operating manuals.