

**Patient Warming System CWS
5000
Commercial Offer**



**Manufacturer: Care Essentials Pty
Ltd (Australia)**

Patient Warming System CWS 5000

The CWS 5000 Patient Warming System is an innovative solution designed to maintain patient normothermia during medical procedures. Its compact and lightweight design makes the CWS 5000 easy to use, while compliance with international safety standards ensures reliable performance in clinical environments.

Features of the CWS 5000 System:

- A dual-speed fan provides two airflow levels, allowing efficient inflation of blankets of various sizes.
- The system has a compact design and lightweight construction without compromising durability.
- User-friendly control interface.
- Status indicators for system filters and integrated temperature sensor system.
- Equipped with retractable hooks for bed mounting and a protective hose cover.
- Includes an external port for calibration and functional testing.

Technical Specifications of the CWS 5000 Patient Warming System:

- Operating voltage: AC 220–240 V
- Power frequency: 50/60 Hz
- Defibrillation protection class (electrical safety): Type BF
- Medical Device Regulation (MDR) classification: Class IIb
- Device temperature settings: 34 / 40 / 43 / 46 °C

Temperature Accuracy and Airflow:

- Temperature accuracy: ± 2 °C (with protective hose cover installed)
- High airflow speed mode: Available via the “HIGH SPEED” button
- (average air temperature measured at hose outlet)

Physical Parameters:

- Dimensions (W×D×H): 290 × 220 × 400 mm
- Device weight: 6 kg

Filtration System:

- Filter type: Bacterial HEPA filter
- Filtration efficiency: Captures particles of 0.2 microns and larger
- Filter service life: 2000 hours



Patient Warming System CWS 5000

*The price does not include additional accessories, blankets, or trolleys.


Prices for accessories are listed below and are calculated based on the packaging quantities of consumables




- **Model: CLM 0101**

Characteristics:


- Full body blanket
- Dimensions (W × L): 100 × 195 cm
- Hose opening: 1
- Delivery quantity: 30 pcs






Supply Specification – Model CLM 0102

- Model: CLM 0102
- Product Type: Lower body blanket
- Dimensions (W × L): 100 × 120 cm
- Hose Opening: 1
- Delivery Quantity: 30 pcs




Supply Specification – Model CLM 0103

- Model: CLM 0103
- Product Type: Upper body blanket (with shoulder coverage)
- Dimensions (W × L): 96 × 130 cm
- Hose Openings: 2
- Delivery Quantity: 30 pcs



Supply Specification – Model CLM 0104

- Model: CLM 0104
- Product Type: Upper body blanket (arms outside)
- Dimensions (W × L): 57 × 196 cm
- Hose Openings: 2
- Delivery Quantity: 30 pcs




Supply Specification – Model CLM 0105

- Model: CLM 0105
- Product Type: Paediatric & neonatal blanket
- Dimensions (W × L): 89 × 95 cm
- Hose Openings: 2
- Delivery Quantity: 10 pcs


Supply Specification – Model CLM 0106

- Model: CLM 0106
- Product Type: Paediatric tubular warming blanket
- Dimensions (W × L): 62 × 125 cm
- Hose Openings: 1
- Delivery Quantity: 10 pcs



Supply Specification – Model CLM 0107

- Model: CLM 0107
- Product Type: Paediatric warming blanket
- Dimensions (W × L): 100 × 135 cm
- Hose Openings: 2
- Delivery Quantity: 10 pcs




Supply Specification – Model CLM 0108

- Model: CLM 0108
- Product Type: Cardiac access warming blanket
- Dimensions (W × L): 100 × 195 cm
- Hose Opening: 1
- Delivery Quantity: 10 pcs



Supply Specification – Model CLM 0109

- Model: CLM 0109
- Product Type: Paediatric full body warming blanket
- Dimensions (W × L): 100 × 120 cm
- Hose Opening: 1
- Delivery Quantity: 10 pcs




Supply Specification – Model CLM 0112




- Model: CLM 0112
- Product Type: Multipurpose air hose sock
- Dimensions: 100 × 14 cm
- Hose Opening: 1
- Delivery Quantity: 30 pcs

Supply Specification – Model CLM 0115



- Model: CLM 0115
- Product Type: Large paediatric underbody warming blanket
- Dimensions (W × L): 100 × 120 cm
- Hose Connection Port: 1
- Quantity per Delivery: 10 units

Supply Specification – Model CLM 0116




- Model: CLM 0116
- Product Type: Small paediatric underbody warming blanket
- Dimensions (W × L): 60 × 90 cm
- Hose Connection Port: 1
- Quantity per Delivery: 10 units

Supply Specification – Model CLM 0117




- Model: CLM 0117
- Product Type: Paediatric lower body warming blanket
- Dimensions (W × L): 60 × 90 cm
- Hose Connection Port: 1
- Quantity per Delivery: 10 units

Supply Specification – Model CLM 0118



- Model: CLM 0118
- Product Type: Abdominal access warming blanket
- Dimensions (W × L): 184 × 195 cm
- Hose Connection Ports: 3
- Quantity per Delivery: 10 units

Supply Specification – Model CLM 0121




- Model: CLM 0121
- Product Type: Multipurpose (Multi-access) warming blanket
- Dimensions (W × L): 100 × 195 cm
- Hose Connection Ports: 1
- Quantity per Delivery: 10 units

Supply Specification – Model CLM 0122






- Model: CLM 0122
- Product Type: Adult tube warming blanket
- Dimensions (W × L): 90 × 217 cm
- Hose Connection Ports: 1
- Quantity per Delivery: 10 units

Supply Specification – Model CLM 0124



- Model: CLM 0124
- Product Type: Surgical access warming blanket (with arm access)
- Dimensions (W × L): 100 × 210 cm
- Hose Connection Ports: 3



Supply Specification – Model CLM 0125

- Model: CLM 0125
- Product Type: Warming blanket for lithotomy & laparoscopy procedures
- Dimensions (W × L): 185 × 295 cm
- Hose Connection Ports: 2
- Quantity per Supply Batch: 10 units

Product Specification – CLM 0127

- Model: CLM 0127
- Product type: Multifunctional blanket (adult)
- Dimensions: 90 × 206 cm
- Number of hose openings: 4
- Quantity for delivery: 10 pcs

Product Specification – CLM 0133

- Model: CLM 0133
- Product type: Under body spinal blanket
- Dimensions: 87 × 198 cm
- Number of hose openings: 1
- Quantity for delivery: 10 pcs

Product Specification – Cocoon CWS 5000

- Model: Cocoon CWS 5000
- Product type: Patient warming system
- Quantity for delivery: 1 pc

Product Specification – CWS 107

- Model: CWS 107
- Product type: Trolley without basket
- Quantity for delivery: 1 pc

Product Specification – CWS 109

- Model: CWS 109
- Product type: Stand with basket
- Quantity for delivery: 1 pc

