

COMMERCIAL PROPOSAL

Medical Beds – TISCOMED

TISCOMED (Turkey)

Model: TBDE 4200

- 4 electric motors
- Patient control panel and nurse control panel integrated into the side rails
- Decentralized brake system (individual wheel locking)



TBDE 4100

- 4 electric motors
- Patient control panel
- Decentralized braking system (individual locking)

+37062812800

medical.trest

TISCOMED TBDE 4100

TECHNICAL SPECIFICATIONS

- 4 electric motors
- Patient control panel
- Backrest, leg section, height adjustment, and Trendelenburg position controlled via patient panel
- Trendelenburg movement is assisted by a gas spring
- 360° rotating wheels with cross brake system
- Brake system with individual pedal control on each wheel
- Mechanical CPR function on back section
- Protective bumpers on all four corners
- Height-adjustable IV pole
- Mattress platform made of ABS plastic
- Urine bag holders on both sides of the bed
- 4-section side rails made of ABS
- Headboard and footboard made of ABS
- Removable headboard and footboard
- CNC-cut mattress with density of 28
- All metal parts coated with electrostatic powder paint for corrosion protection



DIMENSIONS

- Overall dimensions: 1020 × 2200 mm
- Mattress platform: 850 × 1950 mm
- Minimum height: 400 mm
- Maximum height: 800 mm
- Caster diameter: 125 mm



TISCOMED TBDE 4200

TECHNICAL SPECIFICATIONS

- 4 electric motors
- Patient control panel and nurse control panel integrated into the bedside railing
- Adjustment of backrest, leg section, height, and Trendelenburg position via both patient and nurse control panels
- Patient control panel functions can be locked from the nurse control panel
- Fowler and vascular positions available
- Bed features a 360° rotating wheel system with cross-locking
- Decentralized braking system (individual wheel locks)
- Mechanical CPR function on the backrest section
- Protective corner bumpers on all four corners
- Height-adjustable IV pole
- Lying surface (mattress platform) made of ABS plastic
- Urine bag holders on both sides of the bed
- 4-section ABS side rails
- Head and foot boards made of ABS
- Removable head and foot boards
- CNC-cut mattress, density 28
- All metal parts are coated with electrostatic powder paint for corrosion protection

MEDICAL TREST

DIMENSIONS

- Overall dimensions: 1020 × 2200 mm
- Mattress platform size: 850 × 1950 mm
- Minimum height: 400 mm
- Maximum height: 800 mm
- Wheel diameter: 125 mm





**Thank you for your
attention.**

 +37062812800

 medical.trest